



Product Code . JL-E-10514

Zinc Plate with 4mm Plug for Electroscope

Description

If the electroscope is charged and the attached zinc plate is illuminated with UV light from a source such as Mercury lamp, the electroscope will become discharged.

This mounted zinc plate has been designed for use with the electroscope (but can also be used with any electroscope with a 4 mm plug socket) for demonstrating the photoelectric effect.

Jain Scientific Suppliers,
Hargolal Road, Ambala Cantt, Haryana India
Direct Contact Details  +91-8569909696  Jaincolab@gmail.com
 www.jaincolab.in