



Product Code . JL-M-10565

AC Magnetic Field Coil

Description

When AC is passed through the coil, an alternating magnetic field is induced, which can be detected with a search coil.

Jaincolab comprising a 10 turn coil of diameter approximately 15 cm, mounted on a base with 4 mm connecting sockets.

Requires a power supply which can deliver a high current at low voltage.

The current must not be allowed to exceed 6 A.

Jain Scientific Suppliers,
Hargolal Road, Ambala Cantt, Haryana India
Direct Contact Details  +91-8569909696  Jaincolab@gmail.com
 www.jaincolab.in