

Description

Cat. No. : "JSS-527/067"

DESCRIPTION :-

Set of 12 cylindrical bars, of different length & material.

The set includes copper, brass, polyoxymethylene-acetyl, rubber, acrylic, aluminum, polyformaldehyde-Tekaform, PTFE, PVC, oak wood, maple wood, and glass.

The rods are 12mm in diameter and range in length from 37mm to 60mm.

Using this set, your students can explore mass, density, and volume.

Measure the mass using a scale, calculate the volume using a graduated cylinder filled with water, and then calculate the density by dividing the mass by the volume.


Comes in an attractive wooden holder.

Jain Laboratory Instruments Pvt. Ltd,
Hargolal Road, Ambala Cantt, Haryana
India

Direct Contact Details 

+91-8569909696 

sales@jlabexport.com

 www.jlabexport.com